



Measuring mode: positive <input type="checkbox"/> negative <input type="checkbox"/>	Structural formula:
Sample: _____	
Molecular formula: _____	
Molar mass: _____	
Accurate concentration*: _____	
Solvent: _____	
Notes: _____	

Measuring mode: positive <input type="checkbox"/> negative <input type="checkbox"/>	Structural formula:
Sample: _____	
Molecular formula: _____	
Molar mass: _____	
Accurate concentration*: _____	
Solvent: _____	
Notes: _____	

Measuring mode: positive <input type="checkbox"/> negative <input type="checkbox"/>	Structural formula:
Sample: _____	
Molecular formula: _____	
Molar mass: _____	
Accurate concentration*: _____	
Solvent: _____	
Notes: _____	

*Please submit your samples as solution (0.1 – 0.5 mg/mL; at least 100 μ L) in H₂O, MeOH, (iso)propanol, acetonitrile, CH₂Cl₂ (max. 25 % in MeOH) or mixtures from those.